

Instructions for Use

Life Science Kits & Assays



innuPREP Food DNA Kit - PP Mini



Innuscreen
innovative
Sensor
Technology

Order No.:

845-PS-0140016

16 reactions

845-PS-0140096

96 reactions

Publication No.: HB_0140_e_260627

This documentation describes the state at the time of publishing.
It needs not necessarily agree with future versions. Subject to change!

Print-out and further use permitted with indication of source.

© Copyright 2026, IST Innuscreen GmbH

Manufacturer and Distributor:

IST Innuscreen GmbH

Phone +49 30 9489 3380

Robert-Rössle-Straße 10

Fax +49 30 9489 3381

13125 Berlin · Germany

Made in Germany!

info.innu@ist-ag.com

Contents

1	Introduction.....	2
1.1	Intended use.....	2
1.2	Notes on the use of this manual and the kit.....	3
2	Safety precautions.....	4
3	Storage conditions	5
4	Functional testing and technical assistance.....	5
5	Product use and warranty.....	6
6	Kit components	7
6.1	Components included in the kit.....	7
6.2	Components not included in the kit	7
7	Product specifications	7
8	Initial steps before starting.....	8
9	Sample Preparation	9
10	Automated extraction using PurePrep Mini	11
10.1	Prefilling of the DW Plate or the DW Strips	11
10.2	Loading filled Deep Well Plate/Strips to the PurePrep Mini and plug in the Tip Combs	12
11	Troubleshooting.....	14

1 Introduction

1.1 Intended use

The innuPREP Food Kit DNA Kit- PP Mini has been designed for automated isolation of DNA from food samples using the PurePrep Mini instrument. The extraction procedure is based on a new patented chemistry.

The extraction procedure takes place on the magnetic particle processor PurePrep Mini and allows the parallel and flexible extraction of 1 up to 16 samples. The procedure begins with an external lysis step of food samples. The lysed samples are then transferred into the DW Strip or DW Plate (available separately). The following extraction process runs automatically on the PurePrep Mini.

The extraction is based on binding of DNA to surface- modified magnetic particles. After several washing steps the nucleic acids are eluted from the magnetic particles with RNase-free water and are ready to be used in downstream applications. The extraction chemistry in combination with the PurePrep Mini is optimized to get maximum yield and quality.

The kit is intended for use by professional users. The kit has been designed to be used for a wide range of different downstream applications, like amplification reactions and further analytical procedures.











CONSULT INSTRUCTION FOR USE

This package insert must be read carefully before use. Package insert instructions must be followed accordingly. Reliability of results cannot be guaranteed if there are any deviations from the instructions in this package insert.

1.2 Notes on the use of this manual and the kit

For easy reference and orientation, the manual and labels use the following warning and information symbols as well as the shown methodology:

Symbol	Information
	REF Catalogue number.
	Content Contains sufficient reagents for <N> reactions.
	Storage conditions Store at room temperature or shown conditions respectively.
	Consult instructions for use This information must be observed to avoid improper use of the kit and the kit components.
	Expiry date
	Lot number The number of the kit charge.
	Manufactured by Contact information of manufacturer.
	For single use only Do not use components for a second time.
	Note / Attention Observe the notes marked in this way to ensure correct function of the kit and to avoid operating errors for obtaining correct results.

The following systematic approach is introduced in the manual:

- The chapters and figures are numbered consecutively.
- A cross reference is indicated with an arrow (e.g. → „Notes on the use of this manual and the kit“ p. 3).
- Working steps are numbered.

2 Safety precautions

NOTE

Read through this chapter carefully before use to guarantee your own safety and a trouble-free operation.

Follow all the safety instructions explained in the manual, as well as all messages and information, which are shown.

All due care and attention should be exercised in handling the materials and reagents contained in the kit. Always wear gloves while handling these reagents and avoid any skin contact! In case of contact, flush eyes or skin with a large amount of water immediately.



FOR SINGLE USE ONLY!

This kit is made for single use only!

ATTENTION!

Don't eat or drink components of the kit!

The kit shall only be handled by educated personnel in a laboratory environment!

If the buffer bottles are damaged or leaking, wear gloves and protective goggles when discarding the bottles in order to avoid any injuries. This kit can be used with potentially infectious samples. Therefore, all liquid waste must be considered as potentially infectious and must be handled and discarded according to local safety regulation.

Please observe the federal, state and local safety and environmental regulations. Follow the usual precautions for applications using extracted nucleic acids. All materials and reagents used for DNA or RNA isolation should be free of DNases or RNases.

ATTENTION!

Do not add bleach or acidic components to the waste after sample preparation!

NOTE

Emergency medical information in English and German can be obtained 24 hours a day from:

Poison Information Center, Freiburg / Germany

Phone: +49 (0)761 19 240.

For more information on GHS classification and the Safety Data Sheet (SDS) please contact sds.innu@ist-ag.com.

3 Storage conditions

All kit components are shipped at ambient temperature.

Upon arrival, store lyophilized and dissolved **Proteinase K and MAG Suspension** at 4 °C to 8 °C.

All other components of the kit should be stored dry at room temperature (15 °C to 30 °C). When stored at room temperature, the kit is stable until the expiration date printed on the label on the kit box.

If there are any precipitates within the provided solutions solve these precipitates by careful warming. Before every use make sure that all components have room temperature.

4 Functional testing and technical assistance

The IST Innuscreen GmbH guarantees the correct function of the kit for applications as described in the manual. This product has been produced and tested in an ISO 13485 certified facility.

We reserve the right to change or modify our products to enhance their performance and design. If you have any questions or problems regarding any aspects of the **innuPREP Food DNA Kit PP-Mini** or other IST Innuscreen GmbH products, please do not hesitate to contact us. For technical support or further information in Germany please contact info.innu@ist-ag.com. For other countries please contact your local distributor.

5 Product use and warranty

The kit is not designed for the usage of other starting materials or other amounts of starting materials than those referred to in the manual (→ "Product specifications" p. 7). Since the performance characteristics of IST Innuscreen GmbH kits have only been validated for the application described above, the user is responsible for the validation of the performance of IST Innuscreen GmbH kits using other protocols than those described below. IST Innuscreen GmbH kits may be used in clinical diagnostic laboratory systems after the laboratory has validated the complete diagnostic system as required by CLIA' 88 regulations in the U.S. or equivalents in other countries.

All products sold by IST Innuscreen GmbH are subjected to extensive quality control procedures and are warranted to perform as described when used correctly. Any problems should be reported immediately.

NOTE

This kit is for research use only!

6 Kit components

6.1 Components included in the kit

	Σ 16	Σ 96
REF	845-PS-0140016	845-PS-0140096
Lysis Solution CBV	25 ml	150 ml
Proteinase K	2 x for 0.3 ml working solution	2 x for 1.5 ml working solution
MAG Suspension F	0.25 ml	1.1 ml
Binding Solution SBS	8 ml	45 ml
Washing Solution E	12 ml	100 ml
Washing Solution B2 (conc.)	16 ml	80 ml
RNase-free Water	5 ml	30 ml
Manual	1	1

6.2 Components not included in the kit

- 1.5 ml tubes
- 96 %–99.8 % ethanol (molecular biology grade, undenatured)
- ddH₂O; ultrapure for dissolving Proteinase K
- DW Strip / DW Plate / DW Tip Comb (compatible with the PP Mini device)
- Optionally: RNase A solution (10 mg/ml)

7 Product specifications

1. Starting material:
 - Food sample (max. 200 mg)
2. Time for automated extraction protocol on PurePrep Mini:
 - Approx. 15 minutes (excluding external lysis)

8 Initial steps before starting

- Add the indicated volume of ddH₂O to each vial of **Proteinase K**, mix thoroughly and store as described above.

845-PS-0140016	Add 0.3 ml ddH ₂ O to lyophilized Proteinase K.
----------------	--

845-PS-0140096	Add 1.5 ml ddH ₂ O to lyophilized Proteinase K.
----------------	--

- Add the indicated volume of absolute ethanol to **Washing Solution B2 (conc.)** and mix thoroughly. Always keep the bottle firmly closed!

845-PS-0140016	Add 24 ml ethanol to 16 ml Washing Solution B2 .
----------------	---

845-PS-0140096	Add 120 ml ethanol to 80 ml Washing Solution B2 .
----------------	--

- Avoid freezing and thawing of starting material.
- Centrifugation steps should be carried out at room temperature.
- Pre-heat thermal mixer or water bath to 65°C.

9 Sample Preparation

1. Weigh up to 200 mg of food sample. Cut the sample into small piece or homogenise the sample as much as possible before transferring into a 2.0ml tube.
2. Add the recommended amount of Lysis Solution CBV (see table) and 20 µl Proteinase K to each sample and vortex vigorously for 10 seconds. Incubate at 65°C for approximately 60 minutes.

Food Class	Example	Amount of Lysis Solution CBV to add to each sample
Meat Products	Ham, salami	0.8 ml
Tinned food	Fish, meat, sausage	0.8 ml
Milk products	Cheese, yoghurt, chocolate	0.8 ml
Cereals	Flakes, nachos, waffles, cookies, noodles	1.2 ml
Flours	Wheat flour, baking mixes	1.0 ml
Instant products	Instant soups, mashed potato	1.0 ml

NOTE

We recommend using a shaking platform (thermal mixer / water bath / another rocking platform) for a continuous shaking of sample during lysis. Alternatively, vortex the sample every 10 minutes during the incubation. No shaking will reduce the lysis efficiency.

3. Centrifuge the tube at maximum speed for 5 minutes.
4. Transfer the supernatant to a new 1.5 ml reaction tube. If there is floating material above the sample, pierce this carefully with a pi-

Sample Preparation

ette and carefully remove the sample. Avoid aspiration of floating material and/or sediment.

NOTE

To remove RNA from the sample (if necessary), add 2 µl of RNase A solution (10 mg/ml) to the lysed sample, vortex shortly and incubate for 5 minutes at room temperature

5. Check if the sample volume is at least 400 µl. If it is lower, add Lysis Solution CBV up to 400 µl.
6. Proceed with "Automated extraction using PurePrep Mini" on p.12.

10 Automated extraction using PurePrep Mini

10.1 Prefilling of the DW Plate or the DW Strips

	1	2	3	4	5	6	7	8	9	10	11	12
A	Sample 1	→				Eluate 1	Sample 9	→				Eluate 9
B	Sample 2	→				Eluate 2	Sample 10	→				Eluate 10
C	Sample 3	→				Eluate 3	Sample 11	→				Eluate 11
D	Sample 4	→				Eluate 4	Sample 12	→				Eluate 12
E	Sample 5	→				Eluate 5	Sample 13	→				Eluate 13
F	Sample 6	→				Eluate 6	Sample 14	→				Eluate 14
G	Sample 7	→				Eluate 7	Sample 15	→				Eluate 15
H	Sample 8	→				Eluate 8	Sample 16	→				Eluate 16

Fig. 1: Schematic illustration of DW Plate




Fig. 2: Arrangement of the DW Strips in Tray

Cavity of KF96 DW Plate/Strip	Content
Cavity 1 + 7	400 µl lysed sample + 400 µl Binding Solution SBS + 10 µl MAG Solution F
Cavity 2 + 8	600 µl Washing Solution E
Cavity 3 + 9	600 µl Washing Solution B2
Cavity 4 + 10	600 µl Washing Solution B2
Cavity 5 + 11	600 µl Washing Solution B2
Cavity 6 + 12	100 µl – 200 µl RNase-free Water

The prefilling is carried out from left to right as shown in the illustration, Fig. 1. The DW Strips located in the tray are filled in the same way.

10.2 Loading filled Deep Well Plate/Strips to the PurePrep Mini and plug in the Tip Combs

NOTE

- When using a strip (strips), the strip is inserted into the metal tray. In total, a maximum of 8 strips can be used in one extraction-run.
- When working with strips, only every other tip of the tip comb will be used for extraction.

Left tray side: Tip 1, 3, 5, 7

Right tray side: Tip 2, 4, 6, 8.

- It is recommended to mark the tips used for the extraction so that they are not used more than once.

-
1. Select the protocol
"FOOD" and start the run.

2. After finishing the extraction protocol, the Cavities 6 and 12 contain the isolated RNA.
3. Transfer the DNA into a fresh 1.5 ml Tube.

IMPORTANT NOTE

After finishing the extraction protocol, the last cavity of the Plate/Strip contains the isolated DNA. Store the DNA under adequate conditions. We recommend storing the extracted RNA at $-22\text{ }^{\circ}\text{C}$ to $-18\text{ }^{\circ}\text{C}$.

11 Troubleshooting

Problem / probable cause	Comments and suggestions
Little or no total DNA eluted	
Insufficient lysis of starting material	Ensure to use the required volumes of Proteinase K and Lysis Solution CBV .
Eluate volume too high	Decrease the eluate volume. The suggested eluate volume is 150–200 µl. Please note that lowering the eluate volume will not increase the yield proportionally!
Inadequate extraction	Presence of inhibiting substances in the starting material. Please use the kit only for samples that match the requirements declared in "Product specifications".
RNA contamination	
Too much starting material or no RNA digestion with RNase I	Reduce amount of starting material or make sure that RNA digestion done.
Carryover of mag particles	
Eluate contains carryover of magnetic particles	Place the plate on a magnet or centrifuge the plate at maximum speed for 3 minutes. Pipet the supernatant with DNA into a new plate or Elution vessels

iST Innuscreen GmbH
Robert-Rössle-Str.10
13125 Berlin · Germany

Phone +49 30 9489 3380
Fax +49 30 9489 3381

info.innu@ist-ag.com



Innuscreen

innovative
Sensor
Technology