

Instructions for Use

Life Science Kits & Assays



smart DNA prep (a)

Order No.:

845-ASS-2008016 16 reactions
845-ASP-2008016 16 reactions
845-ASS-2008096 96 reactions
845-ASP-2008096 96 reactions

Publication No.: HB_ASP-2008_e_220106

This documentation describes the state at the time of publishing.
It needs not necessarily agree with future versions. Subject to change!

Print-out and further use permitted with indication of source.

© Copyright 2022, IST Innuscreen GmbH

Manufacturer and Distributor:

IST Innuscreen GmbH Phone +49 30 9489 3380
Robert-Rössle-Straße 10 Fax +49 30 9489 3381
13125 Berlin · Germany
Made in Germany! info.innu@ist-ag.com

Contents

1	Introduction.....	3
1.1	Intended use.....	3
1.2	Notes on the use of this manual and the kit.....	4
2	Safety precautions.....	5
3	Storage conditions	6
4	Functional testing and technical assistance.....	7
5	Product use and warranty.....	7
6	Kit components	8
6.1	Included kit components.....	8
6.2	Components not included in the kit.....	8
6.3	Components needed for preparation of bacteria	9
6.4	Components needed for preparation of yeasts	9
7	Initial steps before starting.....	9
8	Product specifications	9
9	Protocols for isolation of DNA	10
9.1	Protocol 1: Isolation from eukaryotic cells	10
9.2	Protocol 2: Isolation from tissue samples	11
9.3	Protocol 3: Isolation from bacteria cell pellets.....	12
9.3.1	Pre-lysis of resuspended starting material...12	
9.3.1.1	Gram-negative bacteria.....	12
9.3.1.2	Gram-positive bacteria.....	12
9.3.1.3	Staphylococcus.....	13
9.3.2	Proteolytic lysis step.....	14
9.4	Protocol 4: Isolation from yeast cell pellets.....	14
9.4.1	Resuspension of starting material.....	14
9.4.2	Pre-lysis of resuspended starting material	14
9.4.3	Proteolytic lysis step.....	15

10	Preparation of Reagent Plates or Reagent Strips	16
10.1	General filling scheme.....	16
10.2	Unpacking of Reagent Plates / Strips and piercing of sealing foil	16
11	Automated extraction using InnuPure C16 / C16 <i>touch</i>	18
11.1	Sample tray of InnuPure C16 / C16 <i>touch</i>	18
11.2	Preparing sample tray of InnuPure C16 / C16 <i>touch</i>	18
11.3	Handling of SmartExtraction Tips.....	20
11.4	Loading the sample to InnuPure C16 / C16 <i>touch</i>	21
11.5	Starting the InnuPure C16	22
11.6	Starting the InnuPure C16 <i>touch</i>	24
12	Troubleshooting.....	27

1 Introduction

1.1 Intended use

The smart DNA prep (a) kit has been designed for automated isolation of high molecular weight genomic DNA (gDNA) from tissue samples, cultured eukaryotic cells, rodent tails, bacteria and yeast cells. The kit utilizes the new SmartExtraction technology invented by IST Innuscreen GmbH.

The procedure starts with the lysis of the starting material. Following lysis of the samples, lysates are transferred into the Reagent Strips or Reagent Plate of the kit, which are already prefilled with all extraction reagents needed for the automated isolation process using a unique 1 ml filter tip in combination with InnuPure C16 / C16 *touch*.

The extraction process is based on adsorption of the genomic DNA to Smart Modified Surfaces inside the tip. After washing, the genomic DNA is eluted from the Smart Modified Surfaces and is ready for use in subsequent downstream applications.

The whole extraction process just needs simple pipetting up and down. The combination of patented, low-salt DC-Technology with patent-pending Smart Modified Surface is optimized to get a maximum of yield and quality.



CONSULT INSTRUCTION FOR USE

This package insert must be read carefully before use. Package insert instructions must be followed accordingly. Reliability of results cannot be guaranteed if there are any deviations from the instructions in this package insert.

1.2 Notes on the use of this manual and the kit

For easy reference and orientation, the manual and labels use the following warning and information symbols as well as the shown methodology:

Symbol	Information
	REF Catalogue number.
	Content Contains sufficient reagents for <N> tests.
	Storage conditions Store at room temperature, unless otherwise specified.
	Consult instructions for use This information must be observed to avoid improper use of the kit and the kit components.
	Expiry date
	Lot number The number of the kit charge.
	Manufactured by Contact information of manufacturer.
	For single use only Do not use components for a second time.
	Note / Attention Observe the notes marked in this way to ensure correct function of the device and to avoid operating errors for obtaining correct results.

The following systematic approach is introduced in the manual:

- The chapters and figures are numbered consecutively.
- A cross reference is indicated with an arrow (e.g. → "Notes on the use of this manual" p. 4).
- Working steps are numbered.

2 Safety precautions

NOTE

Read through this chapter carefully before to guarantee your own safety and a trouble-free operation.

Follow all the safety instructions explained in the manual, as well as all messages and information, which is shown.

All due care and attention should be exercised in handling the materials and reagents contained in the kit. Always wear gloves while handling these reagents and avoid any skin contact! In case of contact, flush eyes or skin with a large amount of water immediately.



FOR SINGLE USE ONLY!

This kit is made for single use only!

ATTENTION!

Don't eat or drink components of the kit!

The kit shall only be handled by educated personnel in a laboratory environment!

If the buffer bottles are damaged or leaking, wear gloves and protective goggles when discarding the bottles in order to avoid any injuries. This kit is to be used with potential infectious human samples. Therefore, all liquid waste must be considered as potentially infectious and must be handled and discarded according to local safety regulation.

Please observe the federal, state and local safety and environmental regulations. Follow the usual precautions for applications using extracted nucleic acids. All materials and reagents used for DNA should be free of DNases or RNases.

ATTENTION!

Do not add bleach or acidic components to the waste after sample preparation!

NOTE

Emergency medical information in English and German can be obtained 24 hours a day from:

Poison Information Center, Freiburg / Germany

Phone: +49 (0)761 19 240.

For more information on GHS classification and the Safety Data Sheet (SDS) please contact sds.innu@ist-ag.com.

3 Storage conditions

All components of the kit are shipped at ambient temperature.

Store lyophilized and dissolved Proteinase K at 4 °C to 8 °C.

All other components of the smart DNA prep (a) should be stored dry at room temperature (15 °C to 30 °C). When stored at room temperature, the kit is stable until the expiration date printed on the label on the kit box. Before every use, make sure that all components have room temperature. If there are any precipitates within the provided solutions dissolve these precipitates by careful warming.

4 Functional testing and technical assistance

The IST Innuscreen GmbH guarantees the correct function of the kit for applications as described in the manual. This product has been produced and tested in an ISO 13485 certified facility.

We reserve the right to change or modify our products to enhance their performance and design. If you have any questions or problems regarding any aspects of the smart DNA prep (a) kit or other IST Innuscreen GmbH products, please do not hesitate to contact us. For technical support or further information in Germany please dial +49 30 9489 3380. For other countries, please contact your local distributor.

5 Product use and warranty

The kit is not designed for the usage of other starting materials or other amounts of starting materials than those, referred to in the manual (→ "Intended use") (→ "Product specifications"). Since the performance characteristics of IST Innuscreen GmbH kits have just been validated for the application described above, the user is responsible for the validation of the performance of IST Innuscreen GmbH kits using other protocols than those described below. IST Innuscreen GmbH kits may be used in clinical diagnostic laboratory systems after the laboratory has validated the complete diagnostic system as required by CLIA' 88 regulations in the U.S. or equivalent regulations required in other countries.

All products sold by the IST Innuscreen GmbH are subjected to extensive quality control procedures and are warranted to perform as described when used correctly. Any problems should be reported immediately.

NOTE

This kit is for research use only!

6 Kit components

6.1 Included kit components

	 16	 96
REF	845-ASS-2008016 ^a 845-ASP-2008016 ^b	845-ASS-2008096 ^a 845-ASP-2008096 ^b
SmartExtraction Tips	16	96
Proteinase K	for 1 x 1.5 ml working solution	for 3 x 1.5 ml working solution
Lysis Solution CBV	10 ml	2 x 25ml
Binding Optimizer	1 ml	5 ml
Reagent Strip L ^a	16 (pre-filled, sealed)	96 (pre-filled, sealed)
Reagent Plate L ^b	2 (pre-filled, sealed)	12 (pre-filled, sealed)
Filter Tips	1 x 16	1 x 96
Elution Tubes (0.65 ml)	16	2 x 48
Elution Caps (Stripes)	2	12
Manual	1	1

6.2 Components not included in the kit

- 1.5 ml and 2.0 ml tubes
- ddH₂O for dissolving **Proteinase K**
- optional RNase A (10 mg/ml)
- 1 x PBS Buffer (137 mM NaCl, 2.7 mM KCl, 10 mM Na₂HPO₄, 1.8 mM KH₂PO₄)
- Piercing Tool, Set (contains Single Piercer and 8 well Piercer; 845-PTS-0000002, IST Innuscreen GmbH, Jena, Germany)

6.3 Components needed for preparation of bacteria

- Lysozyme (stock solution: 10 mg/ml (400 U/μl))
- Mutanolysin (stock solution: 0.4 U/μl)
- Lysostaphin (stock solution: 0.4 U/μl)
- TE-Buffer

Alternatively:

- innuPREP Bacteria Lysis Booster(IST Innuscreen GmbH; 845-KA-1000050)

6.4 Components needed for preparation of yeasts

- Lyticase (stock solution: 10 U/μl)
- Yeast Digest Buffer (50 mM potassium phosphate, 10 mM DTT pH 7.5)

7 Initial steps before starting

- Add 1.5 ml ddH₂O to each vial of **Proteinase K**, mix thoroughly and store as described above.
- Invert the Reagent Plate / Reagent Strips for 3–4 times and thump it onto a table to collect the prefilled solutions at the bottom of the wells.

8 Product specifications

1. Starting material:

- Eukaryotic cells (1×10^5 – 1×10^7)
- Tissue samples (1 mg–100 mg)
- Rodent tail (0.1 cm–1 cm)
- Bacteria cell pellets (1×10^5 – 1×10^9 cells)
- Yeast cell pellets (1×10^5 – 1×10^9 cells)

2. Typical yield:

Not determined. The yield depends on the type and the amount of the starting material

9 Protocols for isolation of DNA

9.1 Protocol 1: Isolation from eukaryotic cells

1. Collect the cells by centrifugation with parameters adequate for the cell type (e.g. 5 minutes at 2,500 x g) and discard the supernatant.
2. Add 200 µl 1 x PBS to the cell pellet and resuspend the pellet as much as possible by intensive pipetting up and down.
3. Add 200 µl Lysis Solution CBV and the needed amount of Proteinase K (see table below), mix vigorously by pulsed vortexing for 5 seconds.

Number of cells	Proteinase K to be added
$1 \times 10^5 - 1 \times 10^6$	20 µl
$1 - 5 \times 10^6$	40 µl

4. Incubate at 55 °C for 30 minutes under continuous shaking.

IMPORTANT

Do not use more starting material as described in "Product specifications" on p. 9!

NOTE

We recommend to use a shaking platform (thermal mixer, water bath or another rocking platform) for a continuous shaking of the sample. Alternatively vortex the sample every 10 minutes during the incubation. No shaking will reduce the lysis efficiency.

NOTE

To remove RNA from the sample (optional) add 1 µl of RNase A solution (10 mg/ml), vortex shortly and incubate for 10 minutes at room temperature. Be sure, that the RNase A is free of DNase-activity.

5. Proceed with "Preparation of Reagent Plates or Reagent Strips" on p. 16.

9.2 Protocol 2: Isolation from tissue samples

1. Cut the starting material into small pieces and place it into a 1.5 ml reaction tube.
2. Add 400 µl Lysis Solution CBV and 40 µl Proteinase K and mix vigorously by pulsed vortexing for 5 seconds.
3. Incubate at 55 °C in a shaking platform until the sample is lysed. Sample lysis time depends on amount and kind of sample. Lysis should be completed within 0.5–3 hours.

IMPORTANT

Do not use more starting material as described in "Product specifications" on p. 9!

NOTE

To remove RNA from the sample (optional) add 1 µl of RNase A solution (10 mg/ml), vortex shortly and incubate for 10 minutes at room temperature. Be sure, that the RNase A is free of DNase-activity.

4. After lysis, centrifuge the 1.5 ml tube at 10,000 x g (12,000 rpm) for 2 minutes to spin down unlysed material.
 5. Use 400 µl of the supernatant for the automated extraction.
-

IMPORTANT

Depending on the sample used, a thin solid phase on the surface of the centrifuged sample can be observed. Avoid the contamination of the sample with this solid phase (e.g. fat).

6. Proceed with "Preparation of Reagent Plates or Reagent Strips" on p. 16.

9.3 Protocol 3: Isolation from bacteria cell pellets

1. Collect the cells by centrifugation with parameters adequate for the cell type (e.g. 10 minutes at 3,000 x g) and discard the supernatant.
2. Resuspend the bacteria cell pellet in **170 µl TE Buffer**. After resuspension start enzymatic pre-lysis as described below. Requirements for pre-lysis depend on the cell type.

9.3.1 Pre-lysis of resuspended starting material

9.3.1.1 Gram-negative bacteria

Although Gram-negative bacteria do not require a pre-lysis step, using Lysozyme (not included in the kit) can enhance the efficiency of lysis.

Using Lysozyme: stock solution of Lysozyme: 10 mg/ml (400 U/µl)

1. Add **20 µl Lysozyme** to the resuspended cells and incubate at 37 °C for 30 minutes under continuous shaking.
2. Proceed with "Proteolytic lysis step" on p. 14.

9.3.1.2 Gram-positive bacteria

Gram-positive bacteria require a pre-lysis step using Mutanolysin and/or Lysozyme (not included in the kit).

Using Lysozyme: stock solution of Lysozyme: 10 mg/ml (400 U/µl)

1. Add **20 µl Lysozyme** to the resuspended cells and incubate at 37 °C for 30 minutes under continuous shaking.

Using Mutanolysin: stock solution of Mutanolysin: 0.4 U/µl

2. Add **5 µl Mutanolysin** to the resuspended cells and incubate at 37 °C for 30 minutes under continuous shaking.
3. Proceed with "Proteolytic lysis step" on p. 14.

NOTE

Lysozyme and Mutanolysin exert synergistic activity. Using both enzymes together will increase the yield of isolated nucleic acids.

Alternatively: use the innuPREP Bacteria Lysis Booster

The innuPREP Bacteria Lysis Booster Kit has been developed for a highly efficient pre-lysis of bacterial cell walls by generating spheroplasts. This new mixture of different enzymes boosts the lysis of all bacteria in particular hard-to-lyse microorganisms like *Streptococcus*, *Lactobacillus*, *Staphylococcus*, *Bacillus* and *Clostridium*.

1. Prepare the enzyme mix according to the manual of the innuPREP Bacteria Lysis Booster.
2. Add 20 µl of the prepared enzyme mix to the sample and vortex it shortly. Incubate the sample for 30 minutes at 37 °C.
3. Proceed with "Proteolytic lysis step" p. 14.

9.3.1.3 Staphylococcus

For lysis of *Staphylococcus* the enzyme Lysostaphin is recommended (not included in the kit).

Using Lysostaphin: stock solution of Lysostaphin: 0.4 U/µl

1. Add 10 µl Lysostaphin to the resuspended cells and incubate at 37 °C for 30 minutes under continuous shaking.
2. Proceed with "Proteolytic lysis step" p. 14.

Alternatively: use the innuPREP Bacteria Lysis Booster

The innuPREP Bacteria Lysis Booster Kit has been developed for a high efficient pre-lysis of bacterial cell walls by generating spheroplasts. This new mixture of different enzymes boosts the lysis of all bacteria in particular hard-to-lyse microorganisms like *Streptococcus*, *Lactobacillus*, *Staphylococcus*, *Bacillus* and *Clostridium*.

1. Prepare the enzyme mix according to the manual of the innuPREP Bacteria Lysis Booster.
2. Add 20 µl of the prepared enzyme mix to the sample and vortex it shortly. Incubate the sample for 30 minutes at 37 °C.
3. Proceed with "Proteolytic lysis step" on p. 14.

9.3.2 Proteolytic lysis step

1. Add **200 µl Lysis Solution CBV** and **30 µl Proteinase K** to the sample and mix vigorously by pulsed vortexing for 5 seconds.
2. Incubate sample for 30 minutes at 55°C and 550 rpm in a shaking platform.

Lysis time of 30 minutes is often sufficient to get enough DNA. If the sample is not clear after 30 minutes prolong the incubation time until the sample is clear.

NOTE

To remove RNA from the sample (optional) add 1 µl of RNase A solution (10 mg/ml), vortex shortly and incubate for 10 minutes at room temperature. Be sure, that the RNase A is free of DNase-activity.

3. Proceed with "Preparation of Reagent Plates or Reagent Strips" p. 16.

9.4 Protocol 4: Isolation from yeast cell pellets

9.4.1 Resuspension of starting material

1. Collect the cells by centrifugation with parameters adequate for the cell type (e.g. 10 minutes with 3,000 x g) and discard the supernatant.
2. Resuspend the yeast cell pellet in **200 µl Yeast Digest Buffer** (→"Components needed for preparation of yeasts" p. 9).
After resuspension start enzymatic pre-lysis as described below.

9.4.2 Pre-lysis of resuspended starting material

For lysis of yeast cells, the enzyme Lyticase is recommended (not included in the kit).

Using Lyticase: stock solution of Lyticase: 10 U/µl

1. Add **10 µl 10 U/µl Lyticase** (not included in the kit) to the resuspended cells and incubate at 37 °C for 30 minutes under continuous shaking.
2. Proceed with "Proteolytic lysis step" p. 15.

9.4.3 Proteolytic lysis step

1. Add **200 µl Lysis Solution CBV** and **30 µl Proteinase K** to the sample and mix vigorously by pulsed vortexing for 5 seconds.
2. Incubate sample for 30 minutes at 55°C and 550 rpm in a shaking platform.

Lysis time of 30 minutes is often sufficient to get enough DNA. If the sample is not clear after 30 minutes prolong the incubation time until the sample is clear.


NOTE

To remove RNA from the sample (optional) add 1 µl of RNase A solution (10 mg/ml), vortex shortly and incubate for 10 minutes at room temperature. Be sure, that the RNase A is free of DNase-activity.

3. Proceed with "Preparation of Reagent Plates or Reagent Strips" on p. 16.

10 Preparation of Reagent Plates or Reagent Strips

10.1 General filling scheme



Cavity 1:	Empty	Cavity 7:	Washing Solution
Cavity 2:	Empty	Cavity 8:	Empty
Cavity 3:	Empty	Cavity 9:	Elution Buffer
Cavity 4:	Binding Solution	Cavity 10:	Empty
Cavity 5:	Washing Solution	Cavity 11:	Washing Solution
Cavity 6:	Washing Solution	Cavity 12:	Empty

10.2 Unpacking of Reagent Plates / Strips and piercing of sealing foil

NOTE

According to transport regulations Reagent Reservoirs are wrapped into plastic bags only when transported by airplane.

A Unpacking of Reagent Reservoirs



Reagent Reservoirs are optional delivered wrapped into plastic bags for transport protection.


Carefully open the overpack of Reagent Reservoirs by using scissors.

B Piercing of sealing foil

NOTE

Invert the Reagent Plates / Reagent Strips 3–4 times and thump it onto a table to collect the pre-filled solutions at the bottom of the wells.

Before using Reagent Plates or Reagent Strips the sealing foil has to be pierced manually. Always wear gloves while piercing of the foil!



Reagent Plates / Reagent Strips are pre-filled with extraction reagents and are sealed with a foil. This foil has to be pierced manually before use, by using the piercing tools (single piercer or 8fold piercer).


NOTE

Keep the Reagent Plates / Reagent Strips in a horizontal position to avoid spilling of the reagents while piercing of the foil.

IMPORTANT NOTE

Open all cavities (one row per sample)!

Using 8 samples in parallel




Using single samples



11 Automated extraction using InnuPure C16 / C16 touch

11.1 Sample tray of InnuPure C16 / C16 touch



No. 1: SmartExtraction and standard filter tips

No. 2: Elution vessels for purified samples

No. 3: Tip block

No. 4: Pressure pad

No. 5: Sample block for reagent plates or adapter for reagent strips

No. 6: Serrated guide rail (C16 *touch*: non-serrated)

No. 7: Adapter for reagent strips

11.2 Preparing sample tray of InnuPure C16 / C16 touch


NOTE

The needed number of Reagent Strips or Reagent Plates is depending on the number of samples, which have to be processed. Don't use more Reagent Strips as number of samples!

1. Move the InnuPure C16 / C16 *touch* sample tray into the Priming Station and fold the holding-down clamp at the sample tray upwards!
2. Place the Reagent Plate or an adapter for the Reagent Strips into the holder of the sample tray. Using Reagent Plates, the notched corner of the Reagent Plate has to align with the colored dot at the holder. Using adapters and Reagent Strips, the colored dot of the adapter has to align with the colored dot at the holder and Reagent Strips have to be inserted in a way that the "AJ" labels are arranged at the side of the adapter, which is more distant from the tip block.


Reagent Plate

The notched corners of the Reagent Plate must point to the colored dot on the holder.



Reagent Strips

Place the Reagent Strips into the adapter. The long tab marked with the label "AJ" must point to the side of the adapter, which is more distant from the tip block.



CAUTION

Both holders have to be equipped with a Reaction Plate or Reagent Strips. If applicable, use an empty or dummy plate for the respective holder.

3. Fold down the holding-down clamp to prevent the Reagent Plates and Reagent Strips to be pulled out of the holder during the extraction process.

4. For each extracted sample place a SmartExtraction Tip and a filter tip in the smaller drill holes of the tip block (→ "Handling of SmartExtraction Tips" p. 20)

NOTE

Extracted high molecular weight DNA from large sample amounts tends to be very viscous. In order to improve the handling of DNA for downstream applications, which don't require high molecular weight DNA, extraction protocols include a homogenization step reducing the fragment size of extracted DNA. If downstream application requires high molecular weight DNA, no standard filter tips may be put in tip row 2. As a result, the eluate will remain in cavities 12 of the reagent plastics at the end of the protocol. Transfer of the eluate into storage tubes (e.g. Elution Tubes with Elution Caps, 1.5 ml reaction tubes) has to be done manually. In order to avoid a loss of DNA integrity pipet carefully with a wide-bore or cut tip.

5. Place the Elution Tubes into the wider drill hole at the edge of the tip block. Empty sample positions do not need to be filled.


NOTE

Especially with the Reagent Strips, make sure that for every Reaction Strip the tips and the elution vessel are in the corresponding positions in the tip block!

IMPORTANT NOTE


Use Elution Tubes (0.65 ml) with corresponding Elution Caps.

11.3 Handling of SmartExtraction Tips



Checking the SmartExtraction Tips.

Make sure that the Smart Modified Material is collected near the outlet of the SmartExtraction Tip. If necessary, flip the tip by finger or edge of table or invert it a few times. The optimal position of the Smart Modified Material inside the tip is shown in the picture on the left side.



11.4 Loading the sample to InnuPure C16 / C16 touch

NOTE

The following step will be done after sample lysis!

1. Prepare the Reagent Plate or Reagent Strips and sample tray according to chapter 10.
-

NOTE

The needed number of Reagent Strips or Reagent Plates is depending on the number of samples, which have to be processed. Don't use more Reagent Strips as number of samples!

2. Transfer the whole sample into the **first cavity** (cavities which are more distant from the tip block) of Reagent Strips or Reagent Plates.
3. Transfer **40 µl** of **Binding Optimizer** to the lysed sample into the **first cavity** of Reagent Strips or Reagent Plates.

11.5 Starting the InnuPure C16

1. Switch on the InnuPure C16 and wait for the device initialization to complete, which is signaled by a beeping sound.
2. Move the loaded sample tray with the Reagent Strips forward into the adapter on the front of the InnuPure C16. The serrated rails at the side of the sample tray must protrude into the grooves of the adapter. After pressing lightly against the tip block the sample tray is automatically pulled into the device.



IMPORTANT – CAUTION

Risk of crushing

Immediately let go of the sample tray once it is being pulled in. Otherwise, there is a risk of your hand being crushed.

3. Start the extraction protocol:
 - Press [SELECT PROTOCOL] in the starting window.
 - Select the desired extraction protocol
 - “SE_Ext_Lysis_C16_01” (32 minutes) or
 - “SE_Ext_Lysis_Fast_C16_01” (21 minutes) or
 - “SE_Ext_Lysis_Sensitive_C16_01” (67 minutes)
- and press [START].

NOTE

For samples with a low or unknown content of DNA always use „Standard protocol” or “Sensitive protocol” for maximum yield. The “Fast protocol” is only recommended for samples with high DNA content in combination with time-critical applications.

4. Enter elution volume and press [OK].

Amount of starting material	Recommended elution volume
< 1 x 10 ⁶ eukaryotic cells	300 µl
1–5 x 10 ⁶ eukaryotic cells	400–500 µl
Bacterial colonies	min. 200 µl
< 1 x 10 ⁸ bacterial cells	min. 200 µl
> 1 x 10 ⁸ bacterial cells	min. 200 µl
< 5 x 10 ⁸ yeast cells	min. 200 µl
> 5 x 10 ⁸ yeast cells	min. 300 µl

5. If needed, choose log file and enter sample ID's, press [OK] or [CANCEL].

NOTE

It is possible to enter sample ID's and to create a run log file. Find more detailed information how to start an extraction protocol using InnuPure C16 (→ user manual p. 37 "6.3.5 Using the sample setup tool").

6. After completion of the protocol press [NEXT] and the sample tray is then automatically moved out of the device.

NOTE

The chosen protocol is performed by the device and after the protocol is finished, the tray with the purified samples will be moved out after pressing [NEXT] and the message 'Program finished' is shown on the screen of the device!

7. Remove the sample tray from the adapter of the InnuPure C16 and move it back into the priming station.
8. After finishing the extraction protocol, the Elution Tubes (0.65 ml) contain the extracted DNA. Close the lids and store the DNA under proper conditions.

NOTE


Store DNA under adequate conditions. We recommend storing the extracted DNA at -22 °C to -18 °C!

11.6 Starting the InnuPure C16 touch

NOTE

The following instructions describe the necessary steps for the start of the InnuPure C16 *touch*. For further features and data entry (e.g. opening templates, entering sample setups, saving projects) refer to the manual of the InnuPure C16 *touch*.

1. Switch on the InnuPure C16 *touch* and the tablet computer. Wait until the home screen of IPextract is displayed on the tablet screen.




NOTE

Home screen of IPextract

2. Choose [TEMPLATES] → [New Template] → [Kit-based].
3. Enter optional information in the tab "General".
4. Choose the tab "Kit Information" and switch the "Technology" to "SmartExtraction"!
5. Choose your desired kit from "Kit Name"!


Automated extraction using InnuPure C16 / C16 touch



NOTE

“Kit Information” tab

6. Enter optional information in the tab “Kit Information”
7. Choose the tab “Extraction” and choose the desired method for “Ethanol Removal” and “Protocol”
 - “Drying” – Ethanol is removed by evaporation
 - “Rinse” – Ethanol is washed away using a special Washing Solution



NOTE

“Extraction” tab

- | | |
|--|-----------------|
| “External Lysis – Drying – 03” | (32 minutes) or |
| “External Lysis – Fast – Drying – 03” | (21 minutes) or |
| “External Lysis – Sensitive – Drying – 03” | (67 minutes) |
| “External Lysis – Rinse – 03” | (28 minutes) or |
| “External Lysis – FAST - Rinse – 03” | (17 minutes) or |
| “External Lysis – Sensitive – Rinse – 03” | (63 minutes) |

NOTE

For samples with a low or unknown content of DNA always use „Standard protocol“ or “Sensitive protocol” for maximum yield. The “Fast protocol” is only recommended for samples with high DNA content in combination with time-critical applications.


NOTE

For most applications, Ethanol Removal by “Drying” is recommended. If the extracted DNA is conceived for very ethanol-sensitive downstream applications (e.g. Droplet PCR), chose the option “Rinse”. “Rinse” can also be selected for time-sensitive preparations, since the protocol saves approx. 6 minutes, but the yield might be lower.

8. Adjust your desired “Eluate Volume” using the slider or the text field. Recommended elution volumes are listed in the table below.

Amount of starting material	Recommended elution volume
< 1 x 10 ⁶ eukaryotic cells	300 µl
1–5 x 10 ⁶ eukaryotic cells	400–500 µl
Bacterial colonies	min. 200 µl
< 1 x 10 ⁸ bacterial cells	min. 200 µl
> 1 x 10 ⁸ bacterial cells	min. 200 µl
< 5 x 10 ⁸ yeast cells	min. 200 µl
> 5 x 10 ⁸ yeast cells	min. 300 µl

9. Choose the tab “Monitoring” and start the protocol by tapping the start button.



NOTE

“Monitoring” tab

10. Follow the instructions displayed on the tablet screen.

11. Completion of the protocol is indicated by a message on the tablet screen. Follow the instructions on the screen to remove the sample tray from the device
12. The Elution Tubes contain the extracted DNA; close the lids and store the DNA under proper conditions.

NOTE

Store the DNA under adequate conditions. We recommend storing the extracted DNA at -22 °C to -18 °C!

12 Troubleshooting

Problem / probable cause	Comments and suggestions
Low amount of extracted DNA	
Insufficient lysis	Increase lysis time. Reduce amount of starting material.
Smart Modified Material not collected near the tip opening	Flip the Pipette Tip by finger or edge of table or invert the Pipette Tip a few times to collect Granulates at the lower part of pipette tip.
Preparation without Binding Optimizer	It is important to add the Binding Optimizer to the Reagent Plastic as described in chapters for handling of the liquid handling platforms. Pay special attention that Binding Optimizer was added to the lysed sample!
High viscosity extracted DNA	
Insufficient amount of Elution Buffer	Elute the DNA with a higher volume of Elution Buffer.
Degraded or sheared DNA	
Old material insufficient	Old material often contains degraded DNA.
RNA contaminations of extracted DNA	RNase A digestion

IST Innuscreen GmbH
Robert-Rössle-Str.10
13125 Berlin · Germany

Phone +49 30 9489 3380
Fax +49 30 9489 3381

info.innu@ist-ag.com