



Innuscreen

innovative
Sensor
Technology

innuPREP PME cfDNA Kit - PP Mini

The innuPREP PME cfDNA Kit - PP Mini is based on a new technology, called: PME – Polymer Mediated Enrichment. The procedure does not start with sample lysis, like commonly used methods or kits. The first step is capturing of free-circulating DNA with a special polymer. Subsequently the captured free-circulating DNA is dissolved in a special lysis buffer and then the DNA is extracted.

The kit combines a unique patented technology for enrichment of biomolecules with a very efficient extraction of DNA using magnetic particles.

Product Name: innuPREP PME cfDNA Kit - PP Mini

Product details

Low Throughput Device: PurePrep Mini

Extract: DNA

Reactions: 16 or 96 (PP Mini)

Sample type/Starting material: Plasma

Specifications:
Starting material

Plasma, liquor, cell culture supernatant (2.5 - 3.5 ml)

Low Throughput Device

PurePrep Mini

Extract

DNA

Reactions

16 or 96 (PP Mini)

Extraction time

90 min (including all steps)

Average yield

Depending on sample quality/quantity

The online shop

Preis: € 192.00

Inhalt: 16 reactions

Please select packing

 ▼