



Innuscreen

innovative
Sensor
Technology

innuPREP Blood DNA Kit - PP Maxi

The innuPREP Blood DNA Kit – PP Maxi has been designed for the automated isolation of genomic DNA from whole blood on the PurePrep Maxi device. The extraction procedure is based on a new kind of chemistry.

The extraction procedure takes place on the magnetic particle processor PurePrep Maxi and allows the parallel and flexible extraction of 1 up to 96 samples.

The extraction process starts with sample lysis on the PurePrep Maxi followed by binding of DNA on the surface of the magnetic particle. After washing the DNA is eluted in Elution Buffer.

The kit is intended for use by professional users. The kit has been designed to be used for a wide range of different downstream applications, like amplification reactions and further analytical procedures.

Product Name: innuPREP Blood DNA Kit - PP Maxi

Product details

High Throughput Device: PurePrep Maxi

Extract: DNA

Reactions: 96, 960 or 9600 (PurePrep Maxi)

Sample type/Starting material: Blood, fresh or frozen

Specifications:
Starting material

200 µl whole blood (EDTA / Citrate)

Average yield

Depending on sample quality/quantity

Extraction time

Appr. 60 minutes for air drying

Appr. 50 minutes for ethanol removal with Washing Solution ER

The online shop

Preis: € 220.00

Inhalt: 96 reactions

Please select packing

▼