



Innuscreen

innovative
Sensor
Technology

innuPREP FFPE DNA Kit - PP Mini

The innuPREP FFPE DNA Kit – PP Mini has been designed for automated isolation of genomic DNA from formalin-fixed, paraffin-embedded (FFPE) samples using the PurePrep Mini device.

The procedure starts with an external lysis step without the need for a deparaffinization step using toxic and hazardous components like octane or xylene. After the external lysis and incubation step the following extraction process runs automatically on the PurePrep Mini. The extraction process is based on binding of the DNA to surface-modified magnetic particles. After washing steps the nucleic acid is eluted from the magnetic particles with Elution Buffer and is now ready to use. The extraction procedure takes place on the magnetic particle processor PurePrep Mini and allows the parallel and flexible extraction of up to 16 samples.

Product Name: innuPREP FFPE DNA Kit - PP Mini

Product details

Low Throughput Device: PurePrep Mini

Extract: DNA

Reactions: 16 or 96 (PP Mini)

Sample type/Starting material: Paraffin samples

Specifications:
Starting material

FFPE tissue samples

Size: approx. $2 \times 5 \mu\text{m}$; optionally more starting material

Average yield

Depending on sample quality/quantity

Extraction time

Preliminary steps: approx. 3.5 hours

Automated Extraction protocol: approx. 30 minutes

The online shop

Price: € 35.00

Content: 16 reactions

Please select packing

▼