



Innuscreen

innovative  
Sensor  
Technology

# innuPREP Blood DNA Kit - PP Maxi

The innuPREP Blood DNA Kit – PP Maxi has been designed for the automated isolation of genomic DNA from whole blood on the PurePrep Maxi device. The extraction procedure is based on a new kind of chemistry.

The extraction procedure takes place on the magnetic particle processor PurePrep Maxi and allows the parallel and flexible extraction of 1 up to 96 samples.

The extraction process starts with sample lysis on the PurePrep Maxi followed by binding of DNA on the surface of the magnetic particle. After washing the DNA is eluted in Elution Buffer.

The kit is intended for use by professional users. The kit has been designed to be used for a wide range of different downstream applications, like amplification reactions and further analytical procedures.

**Product Name:** innuPREP Blood DNA Kit - PP Maxi

## Product details

**High Throughput Device:** PurePrep Maxi

**Extract:** DNA

**Reactions:** 96, 960 or 9600 (PurePrep Maxi)

**Sample type/Starting material:** Blood, fresh or frozen

**Specifications:**  
**Starting material**

200 µl whole blood (EDTA / Citrate)

**Average yield**

Depending on sample quality/quantity

**Extraction time**

Appr. 60 minutes for air drying

Appr. 50 minutes for ethanol removal with Washing Solution ER

## The online shop

Price: € 220.00

Content: 96 reactions

Please select packing

▼