



Innuscreen

innovative
Sensor
Technology

innuPREP Forensic DNA Kit - IPC16

The innuPREP Forensic DNA Kit - IPC16 is an ideal solution to handle smallest forensic samples in an automated process. Thereby the purification of genomic DNA from up to 16 samples in parallel takes place within InnuPure C16 touch. The kit has already been proven with positive result using a large number of starting materials. Each Reagent Strip is prefilled and ready to use for 1 sample.

Those are blood and traces of blood, hair, hair roots and beard stubble, cigarette butts, chewing gum, traces of sperm, swab samples e. g. from ear, as well as finger prints on a variety of surfaces and DNA on a tooth brush. Based on patented technology and the principle of magnetic particle separation, even highly degraded nucleic acids can be recovered. Pre-filled, sealed reagent plastic allows a fast, uncomplicated preparation of the device and isolation of very pure DNA for direct downstream applications.

As an alternative to the standard, pre-filled, sealed extraction kits for InnuPure C16 touch, corresponding lower-cost non-filled kit variants are also available. These kits also contain all the plastics and reagents required for extraction.

Product Name: innuPREP Forensic DNA Kit - IPC16

Product details

Low Throughput Device: InnuPure C16touch

Extract: DNA

Reactions: 16 or 96 (IPC16 - Strip)

Sample type/Starting material: Blood and traces of blood

Hair

Fingernails

Stamps and envelopes

Cigarette butts

Chewing gum

Swabs

Sperm

Bone meal

Specifications:

Extraction of high quality DNA for immediate downstream applications

For usage of InnuPure C16 touch and up to 16 samples in parallel

Starting material

Blood and traces of blood

Hair, hair roots and beard stubble

Finger nails

Stamps and envelopes

Cigarette butts, chewing gum

Swab samples and fingerprints taken from surfaces, ear swab, tooth brush

Traces of sperm, bone meal

Extraction time

Lysis: approx. 120 minutes (external)

InnuPure C16 touch protocol: approx. 38 - 47 minutes

Average yield

Depends on the sample and the amount used

Average purity

1,8 - 2,0

The online shop

Preis: € 152.25

Inhalt: 16 reactions

Please select packing

16 reactions ▼