



Innuscreen

innovative
Sensor
Technology

innuSPEED Lysis Tube E

Lysis Tube E is a 2.0 mL tube with screw cap (orange) containing optimized ceramic beads (2.4–2.6 mm) for efficient disruption of difficult samples as insects, dried tissues or plants using SpeedMill or other homogenizers.

Product Details

- Fast, efficient homogenization of robust samples for isolating nucleic acids or proteins
- 0.5 mL and 2.0 mL reaction vessels with beads
- Flexibility thanks to a variety of bead materials and sizes
- Ideal for use with the SpeedMill or other commercially available homogenizers

innuSPEED lysis tubes include a variety of 0.5 mL and 2.0 mL reaction vessels with skirted bases and screw caps. Pre-filled with beads in various sizes, materials and quantities, innuSPEED lysis tubes are ideally suited for homogenizing an exceptionally wide variety of starting materials (e.g., plants, tissues, cells, etc.).

As a general rule for choosing the right matrix: the smaller the sample, the smaller the beads.

Product Name: innuSPEED Lysis Tube E

Product details

Lysis Tubes: Ceramic beads 2.4–2.8 mm

The online shop

Preis: € 31.50

Inhalt: 10 reactions

Please select packing

10 reactions ▼