



Innuscreen

innovative
Sensor
Technology

innuSPEED Lysis Tube P

Lysis Tube P is a 0.5 mL tube with screw cap (green) containing optimized ceramic beads (2.4–2.8 mm) for efficient disruption for difficult samples as insects, dried tissues or plants using SpeedMill homogenizers.

Product Details

- Fast, efficient homogenization of robust samples for isolating nucleic acids or proteins
- 0.5 or 2.0 mL reaction vessels with beads
- Flexibility thanks to a variety of bead materials and sizes
- Ideal for use with the SpeedMill or other commercially available homogenizers

innuSPEED lysis tubes include a variety of 0.5 mL and 2.0 mL reaction vessels with skirted bases and screw caps. Pre-filled with beads in various sizes, materials and quantities, innuSPEED lysis tubes are ideally suited for homogenizing an exceptionally wide variety of starting materials (e.g., plants, tissues, cells, etc.).

As a general rule for choosing the right matrix: the smaller the sample, the smaller the beads.

Product Name: innuSPEED Lysis Tube P

Product details

Lysis Tubes: Ceramic beads 2.4–2.8 mm

The online shop

Preis: € 31.50

Inhalt: 10 reactions

Please select packing

10 reactions ▼