



Innuscreen

innovative
Sensor
Technology

50x inNucleotide Mix

IST Innuscreen's ready-to-use nucleotide mix provides highest quality of deoxynucleotides. All dNTPs are ultra pure (> 98 %) and quality checked by a set of PCR, RT-PCR and Klenow reactions. It is supplied as a 50-fold concentrated mix of ultra pure dATP, dCTP, dGTP and dTTP with 12.5 mM each. The total amount of dNTP in each tube is 25 µmol (6.25 µmol of each dNTP).

Specifications

Concentration

12.5 mM pro dNTP

Quality Control

The deoxynucleotide solution has a purity of > 98 % and is functionally tested by real-time amplification of 30 kb DNA fragments.

Packaging Size

2 × 25 µmol Store at -20 °C, avoid frequent thawing and freezing.

Product Name: 50x inNucleotide Mix

The online shop

Preis: € 8.00

Inhalt: 10 reactions

Please select packing

10 reactions ▼