



Innuscreen

innovative
Sensor
Technology

innuPREP Blood DNA Kit - PP Mini

The innuPREP Blood DNA Kit – PP Mini has been designed for the automated isolation of genomic DNA from whole blood using the PurePrep Mini device.

The extraction procedure takes place on the magnetic particle processor PurePrep Mini and allows the parallel and flexible extraction of up to 16 samples. The extraction process starts with sample lysis on the PurePrep Mini followed by binding of DNA on the surface of the magnetic particle. After washing the DNA is eluted in Elution Buffer.

Product Name: innuPREP Blood DNA Kit - PP Mini

Product details

Low Throughput Device: PurePrep Mini

Extract: DNA

Reactions: 16 or 96 (PP Mini)

Sample type/Starting material: Blood, fresh or frozen

Specifications:
Starting material

Whole blood (up to 200µl)

Fresh or frozen whole blood samples

Stabilized with EDTA or Citrate

Average yield

Depending on sample quality/quantity

Extraction time

58 min

The online shop

Preis: € 35.00

Inhalt: 16 reactions

Please select packing

16 reactions

