



Innuscreen

innovative
Sensor
Technology

innuPREP Plate Set - PP Mini

Consumables for the magnetic bead-based innuPREP kits of the PurePrep Mini

16 rxn - 1 DW plate and 2 tip combs

96 rxn - 6 DW plates and 12 tip combs

Product Name: innuPREP Plate Set - PP Mini

Product details

Low Throughput Device: PurePrep Mini

Extract: DNA & RNA

The online shop

Preis: € 13.00

Inhalt: 16 reactions

Please select packing ▼