



Innuscreen

innovative
Sensor
Technology

innuPREP Strip Set - PP Mini

Consumables for the magnetic bead-based innuPREP kits of the PurePrep Mini

16 rxn - 16 DW strips and 4 tip combs

96 rxn - 96 DW strips and 24 tip combs

Product Name: innuPREP Strip Set - PP Mini

Product details

Low Throughput Device: PurePrep Mini

Extract: DNA & RNA

The online shop

Preis: € 22.00

Inhalt: 16 reactions

Please select packing ▼