



Innuscreen

innovative
Sensor
Technology

innuPREP Leaf & Seed Plate Set - PP Maxi

Consumables for the extraction of 1 up to 960 samples in combination with the PurePrep Maxi device and the innuPREP Leaf & Seed (ultra-fast) - PP Maxi Kit.

The extraction procedure takes place on the magnetic particle processor PurePrep Maxi and allows the parallel and flexible extraction.

960 reactions are including: 10 Tip Combs (96er); 10 Deep Well Plates (2 ml); 50 Flat Well Plates

The kit is intended for use by professional users.

Product Name: innuPREP Leaf & Seed Plate Set - PP Maxi

Product details

High Throughput Device: PurePrep Maxi

Extract: DNA & RNA

Reactions: 96, 960 or 9600 (PurePrep Maxi)

The online shop

Preis: € 300.00

Inhalt: 960 reactions