



Innuscreen

innovative
Sensor
Technology

innuSPEED Lysis Tube S

Lysis Tube S is a 2.0 mL tube with screw cap (yellow) containing special small ceramic beads (0.4–0.6 mm) for efficient disruption of stool, bacteria and fungi samples using SpeedMill and other homogenizers.

Product Details

- Fast, efficient homogenization of robust samples for isolating nucleic acids or proteins
- 0.5 or 2.0 mL reaction vessels with beads
- Flexibility thanks to a variety of bead materials and sizes
- Ideal for use with the SpeedMill or other commercially available homogenizers

innuSPEED lysis tubes include a variety of 0.5 mL and 2.0 mL reaction vessels with skirted bases and screw caps. Pre-filled with beads in various sizes, materials and quantities, innuSPEED lysis tubes are ideally suited for homogenizing an exceptionally wide variety of starting materials (e.g., plants, tissues, cells, etc.).

As a general rule for choosing the right matrix: the smaller the sample, the smaller the beads.

Product Name: innuSPEED Lysis Tube S

Product details

Lysis Tubes: Ceramic beads 0.4-0.6 mm

The online shop

Price: € 31.50

Content: 10 reactions

Please select packing

10 reactions ▼