



Innuscreen

innovative  
Sensor  
Technology

# **innuDETECT Salmonella spp. Assay - Food Pathogen**

TaqMan®-based innuDETECT Food Pathogen Assays enable specific analysis for a range of pathogens encountered in food control applications.

Are you looking for a fast and sensitive method for the identification of microorganisms in food samples? Thanks to developments in molecular biological techniques, IST Innuscreen's solutions minimize the analysis time with minimal personnel effort and reduce the costs incurred while taking into account the quality requirements of today's daily laboratory routines. Benefit from the new TaqMan®-based innuDETECT Food Pathogen Assays, which enable specific analysis for a range of pathogens encountered in food control applications. DNA extraction using our novel SmartExtraction technology and subsequent application of IST Innuscreen's innuDETECT Food Pathogen Assays, based on lyophilized master mixes, make the daily use easy and secure. The reproducible real-time PCR results are characterized by very good efficiency, detected in the FAM channel. The artificial internal control integrated in the assays can be used as extraction as well as amplification control and is detected in the HEX channel.

## **Specifications**

- Verified for use after manual and automated nucleic acid extraction
- Lyophilized reagents for easy PCR set-up, to minimize errors, for improved efficiency and for shipping at room temperature
- Internal control included
- Fast analysis due to optimized PCR chemistry: assay time 1 h

## **Detection principle**

Real-Time PCR (TaqMan)

Target-DNA (FAM); Internal Control (HEX)

**Starting material**

Pathogen DNA isolated from pellets following stomacher cultivation

**Detection Time**

1 h

**Sensitivity**

5 Copies/PCR reaction

**Product Name:** innuDETECT Salmonella spp. Assay

**Product details**

**Reactions:** 24 or 96/Food Pathogen Assays

**The online shop**

Price: € 126.00

Content: 24 reactions

Please select packing  ▼