



Innuscreen

innovative  
Sensor  
Technology

# deltaPREP Blood DNA Mini Kit (MDX)

The deltaPREP Blood DNA Mini Kit (MDX) is a highly efficient tool for directly isolating DNA from whole blood samples of up to 400  $\mu$ l.

The optimized kit for extraction of DNA from blood guarantees extremely high-quality DNA that can be used immediately for photometric determinations, PCR or other downstream applications. The deltaPREP Blood DNA Mini Kit (MDX) is very fast and easy to use. Isolation of genomic DNA – free of inhibitors or impurities – can be completed in just 24 minutes.

**Product Name:** deltaPREP Blood DNA Mini Kit (MDX)

## Product details

**Extract:** DNA

**Reactions:** 10 / 50 / 100 / 250 (either or)

**Sample type/Starting material:** Blood, fresh or frozen (RNA only fresh)

**Specifications:**  
**Starting material**

- Whole blood samples (up to 400  $\mu$ l)
- Fresh or frozen blood
- Stabilizers: EDTA or citrate

### **Average yield**

Depends on sample and the used volume

Up to 30  $\mu$ g gDNA

### **Extraction time**

Approx. 24 minutes, including lysis

### **Binding capacity**

Column binding capacity: > 50  $\mu$ g gDNA

## Average purity

1,7 - 2,0

## The online shop

Price: € 42.00

Content: 10 reactions

Please select packing

 ▼