



Innuscreen

innovative
Sensor
Technology

innuSPEED Lysis Tube H

Lysis Tube H is a 0.5 mL tube with screw cap (violet) containing optimized steel beads (8 x 3.5 mm) for efficient disruption of hard samples as wood, seed or rice grain using SpeedMill or other homogenizers.

Product Details

- Fast, efficient homogenization of robust samples for isolating nucleic acids or proteins
- 0.5 or 2.0 mL reaction vessels with beads
- Flexibility thanks to a variety of bead materials and sizes
- Ideal for use with the SpeedMill or other commercially available homogenizers

innuSPEED lysis tubes include a variety of 0.5 mL and 2.0 mL reaction vessels with skirted bases and screw caps. Pre-filled with beads in various sizes, materials and quantities, innuSPEED lysis tubes are ideally suited for homogenizing an exceptionally wide variety of starting materials (e.g., plants, tissues, cells, etc.).

As a general rule for choosing the right matrix: the smaller the sample, the smaller the beads.

Product Name: innuSPEED Lysis Tube H

Product details

Lysis Tubes: Steel beads 3.5 mm

The online shop

Price: € 31.50

Content: 10 reactions

Please select packing

10 reactions ▼