

# innuSPEED Lysis Tube W

**Lysis Tube W** is a 2.0 mL tube with screw cap (red) containing a mixture of ceramic and steel beads (1.4-1.6 mm and 3.5 mm) for efficient disruption of stool, bacteria, fungi and soil samples as well as cell cultures and spores using SpeedMill or other homogenizers.

## Product Details

- Fast, efficient homogenization of robust samples for isolating nucleic acids or proteins
- 0.5 or 2.0 mL reaction vessels with beads
- Flexibility thanks to a variety of bead materials and sizes
- Ideal for use with the SpeedMill or other commercially available homogenizers

innuSPEED lysis tubes include a variety of 0.5 mL and 2.0 mL reaction vessels with skirted bases and screw caps. Pre-filled with beads in various sizes, materials and quantities, innuSPEED lysis tubes are ideally suited for homogenizing an exceptionally wide variety of starting materials (e.g., plants, tissues, cells, etc.).

As a general rule for choosing the right matrix: the smaller the sample, the smaller the beads.

**Product Name:** innuSPEED Lysis Tube W

## Product details

**Lysis Tubes:** Ceramic and steel beads mixed (1.4- 1.6 mm and 3.5 mm)

## The online shop

Price: € 31.50

Content: 10 reactions

Please select packing

 ▼